

Minor Surgery & Diagnostics

Outpatient® II

Outpatient® II is the industry leader for diagnostic and procedure lights. With minimal heat output and high color temperature, the Outpatient® II is the perfect low maintenance, top performing light for nearly every application in the hospital, surgical suite or physician's office.

The Outpatient® II utilizes three powerful 50 watt halogen lamps that reduce dangerous shadows. Its drift-free design and heavy duty construction make it your best choice in medical task lighting.

(Formerly known as Trinity)



FEATURES



Outpatient® II single ceiling



Outpatient® II double ceiling



Outpatient® II floor stand



Outpatient® II wall mount

Specifications

(Specifications subject to change, consult web site for any revisions)

- Housing:** Formed Lexan housing and handle; 16GA welded steel arms and 11 GA welded steel extension arm; durable swivel head mount; die cast zinc yoke.
- Switching:**
 - On/off switching located in the pivot arm for all models.
 - Floor stand: on/off foot switch.
- Lamp:** Single Outpatient® II head includes (3) 50 watt halogen lamps wired in parallel for continuous operation; rated for 2000 hours; 3300°K, 97 CRI; output 8,000 foot candles at 24 inches *(lamps included)*.
- Lamp Shield:** Heat-filters and dichroic reflectors.
- Focusing:** Adjustable focus with rotation of center autoclavable handle.
- Mounting:** Available in single and double head fixed and track & trolley systems, wall mount and floor stand models.
- Reach & Sweep:**
 - Ceiling mount vertical sweep: 28 inches.
 - Ceiling mount horizontal radius: 180°.
 - Wall mount vertical sweep: 32.5 inches to 63 inches from the floor.
 - Wall mount horizontal radius: 180°.
- Finish:** Housing, yoke, arms and handle are baked enamel powder coated light grey.
- Handle:** One autoclavable handle comes standard with every unit; disposable handle covers are available.
- CA OSHPD:** California Office of Statewide Health Planning and Development: Mounting system meets California seismic codes.
- Manufacturer:** Outpatient® II is manufactured by Burton Medical, Inc. Healthcare Lighting® is a distributor for this product.
- Listings:** UL Listed 60601-1/IEC 60601-2-4 certified.



Single Outpatient® II on track & trolley system



Double Outpatient® II on track & trolley system



ORDERING INFORMATION

Example: HMS301-TT-120-6130

Series ⁽¹⁾

		Description.....			
HMS301-C-120	Single Outpatient® II head	(Ceiling Mount-120v)	(3) 50 watt halogen lamps (included)		
HMS301-F-120	Single Outpatient® II head	(Floor Stand Mount-120v)	(3) 50 watt halogen lamps (included)		
HMS301-W-120	Single Outpatient® II head	(Wall Mount-120v)	(3) 50 watt halogen lamps (included)		
HMS301-TT-120	Single Outpatient® II head	(Track/Trolley System-120v)	(3) 50 watt halogen lamps (included)		
HMS302-C-120	Double Outpatient® II head	(Ceiling Mount-120v)	(6) 50 watt halogen lamps (included)		
HMS302-TT-120	Double Outpatient® II head	(Track/Trolley System-120v)	(6) 50 watt halogen lamps (included)		

Options

8100	Box of 25 SteriHandle™ disposable covers				
8465	Replacement autoclavable handle				
6130	Replacement 50 watt, 120v lamps (pack of 3)				

Ordering Notes: (1) Product available in 230v (International Sales Only)

