





# MED-GAS SERVICE CONSOLE



## Specifications

- Construction:** 20GA galvanized steel housing and ends; heavy gauge aluminum face plate.
- Electrical:** Standard duplexes (ivory) and emergency duplexes (red) shall be 20 amp, Hospital grade type.
- Gas Outlets:** All major manufacturers, DSS and Quick Connect type fittings; both provision-only or factory sourcing are available.
- Slide Brackets:** Will accommodate standard vacuum bottles.
- Nurse Calls:** Provisions only; specify brand and model number for proper dimensions.
- Telephone:** Provisions only for phone jack; Finish plate and internal components by others.
- Special:** Contact factory for any special device not listed
- Knock-outs:** Each gang shall have a 7/8 inch knock out to provide access to all devices gangs.
- Blanks:** All blanks consist of an empty gang area; face covers are devoid of all openings, out lines or markings; blanks are used to space out devices to a desired spacing.
- Standard Finish:** Housing back box is galvanized; face plate is brush aluminum.
- Listings:** UL/C-UL listed.



Specifications subject to change. Please visit our web site for updated information.

## Ordering Information

When ordering or requesting a quotation, please submit a console map using the below nomenclature on this sheet or from your own engineering diagrams.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	16	17	18	19
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Common Gas Outlets

- O2 Oxygen
- A Air
- V Vacuum
- NO2 Nitrous Oxide
- N2 Nitrogen
- CO2 Carbon dioxide

### Common Device Accessories

- DX Duplex receptacle, 20a Ivory
- EDX Duplex receptacle, 20a Red
- TS Toggle switch, 20a, Ivory
- ETS Toggle switch, 20a, Red
- SL Vacuum bottle slide bracket
- GRD Hospital grade grounding jack
- TEL Provisions for Telephone jack

### Nurse Call Provisions

- NC-11 1-gang Nurse Call
- NC-2 2-gang Nurse Call
- NC-3 3-gang Nurse Call
- NC-4 4-gang Nurse Call
- CB-1 1-gang Code Blue
- CB-2 2-gang Code Blue

### Required Information

- Nurse call brand: \_\_\_\_\_
- Nurse call MDL #: \_\_\_\_\_
- Code blue brand: \_\_\_\_\_
- Code blue Mdl #: \_\_\_\_\_
- Gas outlet type: \_\_\_\_\_
- Gas outlet brand: \_\_\_\_\_

### ATTENTION

All Console orders must have approved, factory-engineered drawings before release. Lead-times may vary. Contact the factory for delivery details.